


# AED Automated External Defibrillator

HEART  GUARDIAN



Technology for Life



## HR-503T (Trainer) HR-503 (AED)

**Radianqbio AED (Automated External Defibrillator)** helps recover normal cardiac function upon acute heart disease, such as acute cardiac arrest, by providing cardiac function analysis and, if necessary, electroshock.

Heart Guardian™ as an AED, provides precise analysis of victims' body conditions and configures exact amount of shocks needed, promising safety and quicker responses.





	HR-503	HR-503T
Product Type	AED Life Mode	AED Trainer
ECG Analysis	8.5 Sec.	
Hi- Cap Charging	8 Sec.	
Shock Output	150J(Adults), 50J(Kids)	
Data Storage	40 Min. (Embedded memory)	
Data Transfer Protocol	Bluetooth 2.0	
Control Buttons	Power, Shock, Language, Adult/Child	Power, Shock, Language, Adult/Child
Battery Type	Li-Mn(12V DC 4.2A)	Ni-MH(9.6V DC 2A)
Shock Frequency	200, at least	
Pad Type	Pre-connectyde Electrode Pads, disposable and adhesive	Pads for training purpose (Washable Pads)
Dimension	224(W) x 312(L) x 92(H)mm	224(W) x 312(L) x 92(H)mm
Weight	2.3 kg	1 kg
Warranty	5 Yrs. (2 Yrs. for pads/ 4 Yrs. for Batteries)	1 Yrs. (6 months for Batteries)
Paid Repair Service	10 Yrs	10 Yrs