

TREADMILLS

Heavy Duty Medical Treadmill for Exercise Testing

THE BEST CHOICE FOR
MEDICAL USE



TM8

TM9

NORAV
medical
ECG Excellence

Features

Triple thickness designed to guarantee maximum energy absorption making the use as comfortable as possible.

Non-slip rubber coating for maximum patient safety, contributing at the same time to the structural and mechanical rigidity of the system.

Dimensions

Length
Width
Height
(with treadmill at zero
inclination)

TM8

205 cm
80 cm
112 cm

TM9

230 cm
89 cm
115 cm

		TM8		TM9	
Running surface		50 x 150cm		58 x 170cm	
Distance of the running surface from the floor		17cm		17cm 93cm	
Distance of the handrails from the running surface		93cm			
Weight		165kg		220kg	
Patient capacity		200kg		200kg	
Speed range		0 - 20 km/h		0 - 25 km/h	
Speed increments		0,1 km/h		0,1 km/h	
Incline range		0 - 25%		0 - 25%	
Incline increments		0,5 %		0,5 %	
Speed motor		Asynchronous three phase 2,0HP, AC (nominal power)		Asynchronous three phase 3,0HP, AC (nominal power)	
Elevation motor		90 Watt DC (nominal power)		90 Watt DC (nominal power)	
Running belt		Impact absorbing (Thickness 3,1mm)		Impact absorbing (Thickness 3,1mm)	
Colors					
Surfaces		embossed white		embossed white	
Handrails		black		black	
Belt		black		black	
Safety devices		2 red stop switches		2 red stop switches	
		Safety lanyard (optional)		Safety lanyard (optional)	
		Pediatric handrails (optional)		Pediatric handrails (optional)	
Data Input/Output					
Serial port		RS232 asynchrone		RS232 asynchrone	
Connector		DB9 F		DB9 F	
Baud rate		4800 baud rate		4800 baud rate	
Operation		Remote device interface		Remote device interface	
Utilisation conditions					
Temperature		from 10 to 40°		from 10 to 40°	
Humidity		<85%		<85%	
Storage and transport condition					
Temperature		from -10 to 60°		from -10 to 60°	
Humidity		<90%		<90%	
Packaging		Carton box bound with wooden pallet		Carton box bound with wooden pallet	
		218x92x68 h cm / Gross weight kg 200		244x100x73 h cm / Gross weight Kg 260	
Technical specifications					
Power supply		Standard		Standard	
Voltage		230V 50/60Hz		230V 50/60Hz	
Power		1,5KVA		2,2KVA	
Power supply cord		2P+T, 3 x 2,5mm ² , lenght 3,0m		2P+T, 3 x 2,5mm ² , lenght 3,0m	
					



* With standard handrails



* Optional addition
Control Panel



* Optional pediatric handrails



* With optional long handrails

Optional Accessories

- Pediatric handrails
- Safety strain switch (auto-stop function)

Touch Control Panel for use independent of PC

NORAV
medical
ECG Excellence

www.noravmedical.com | For additional information please contact us at info@norav.com